



VWF Antibody, Biotin conjugated

Product Code	CSB-PA12779D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04275
Immunogen	Recombinant Human von Willebrand factor protein (1-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	von Willebrand factor (vWF) [Cleaved into: von Willebrand antigen 2 (von Willebrand antigen II)], VWF, F8VWF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VWF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.