



# Riboflavin-binding protein Antibody

<b>Product Code</b>	CSB-PA360879YA01CH
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02752
<b>Immunogen</b>	Recombinant Chicken Riboflavin-binding protein (18-225AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Chicken
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Riboflavin-binding protein (RBP) [Cleaved into: Riboflavin-binding protein, plasma form; Riboflavin-binding protein, yolk major form; Riboflavin-binding protein, yolk minor form]
<b>Immunogen Species</b>	Chicken
<b>Research Area</b>	Others