



gag Antibody, Biotin conjugated

Product Code	CSB-PA356017DD01MHK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03332
Immunogen	Recombinant Moloney murine leukemia virus Gag polyprotein protein (216-478AA)
Raised In	Rabbit
Species Reactivity	Moloney murine leukemia virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gag polyprotein (Pr65gag) (Core polyprotein) [Cleaved into: Matrix protein p15 (MA); RNA-binding phosphoprotein p12 (pp12); Capsid protein p30 (CA); Nucleocapsid protein p10-Gag (NC-gag)], gag
Immunogen Species	Moloney murine leukemia virus
Research Area	Others
Target Names	gag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.