



FTL Antibody, HRP conjugated

Product Code	CSB-PA11407B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02792
Immunogen	Recombinant Human Ferritin light chain protein (2-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ferritin light chain (Ferritin L subunit), FTL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FTL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.