




AFP Antibody

Product Code	CSB-PA09989A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02771
Immunogen	Recombinant Human Alpha-fetoprotein protein (19-609AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), AFP, HPAFP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	AFP

Image



Western Blot

Positive WB detected in: 293T whole cell lysate
 All lanes: AFP antibody at 3.2µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 69 kDa
 Observed band size: 69 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.