




TTR Antibody

Product Code	CSB-PA12077A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02766
Immunogen	Recombinant Human Transthyretin protein (21-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TTR

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA12077A0Rb at dilution of 1:10