



HBG1 Antibody, Biotin conjugated

Product Code	CSB-PA010155HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69891
Immunogen	Recombinant Human Hemoglobin subunit gamma-1 protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit gamma-1 (Gamma-1-globin) (Hb F Agamma) (Hemoglobin gamma-1 chain) (Hemoglobin gamma-A chain), HBG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HBG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.