



mpt64 Antibody, FITC conjugated

Product Code	CSB-PA358713YC01FSG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WIN9
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunogenic protein MPT64 (Antigen MPT64), mpt64
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Research Area	Others
Target Names	mpt64
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.