



SFTPC Antibody

Product Code	CSB-PA021174LA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15783
Immunogen	Recombinant Bovine Pulmonary surfactant-associated protein C protein (25-58AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pulmonary surfactant-associated protein C (SP-C) (Pulmonary surfactant-associated proteolipid SPL(Val)), SFTPC, SFTP2
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	SFTPC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.