




ADIPOQ Antibody

Product Code	CSB-PA07957A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15848
Immunogen	Recombinant Human Adiponectin protein (19-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipose most abundant gene transcript 1 protein) (apM-1) (Gelatin-binding protein), ADIPOQ, ACDC ACRP30 APM1 GBP28
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ADIPOQ

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA07957A0Rb at dilution of 1:100