



ESR1 Antibody, Biotin conjugated

Product Code	CSB-PA007830YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03372
Immunogen	Recombinant Human Estrogen receptor protein (10-591AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1), ESR1, ESR NR3A1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ESR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.