



VP3 Antibody, FITC conjugated

Product Code	CSB-PA860325DC01CID
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99152
Immunogen	Recombinant Chicken anemia virus Apoptin protein (1-121AA)
Raised In	Rabbit
Species Reactivity	Chicken anemia virus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptin, VP3
Immunogen Species	Chicken anemia virus
Research Area	Others
Target Names	VP3