



ALB Antibody, FITC conjugated

Product Code	CSB-PA00950G1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02769
Immunogen	Serum Albumin purified from Bovine serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (BSA) (allergen Bos d 6), ALB
Immunogen Species	Bos taurus (Bovine)
Research Area	Cardiovascular
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.