




# Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] (poxB)

<b>Product Code</b>	CSB-EP356996ENV
<b>Abbreviation</b>	Recombinant E.coli poxB protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P07003
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MKQTVAAAYIAKTLESAGVKRIWGVTDGDSLNGLSDSLNRMGTIIEWMSTRHEEV AAFAAGAEALQSGELAVCAGSCGPGNLHLINGLFDCHRNHVPVLAIAAHIPSSE IGSGYFQETHPQELFRECSHYCELVSSPEQIPQVLAIAMRKAVLNRGVSVVLP GDVALKPAPEGATMHWYHAPQPVVTPEEEELRKLALQLLRYSSNIALMCGSGC AGAHKELVEFAGKIKAPIVHALRGKEHVEYDNPYDVGMTGLIGFSSGFHTMMN ADTLVLLGTQFPYRAFYPYPTDAKIIQIDINPASIGAHSKVDMALVGDIAKSTLRALLP LVEEKADRKFLDKALEDYRDARKGLDDLAKPSEKAIHPQYLAQQISHFAADDAI FTCDVGTPTVWAARYLKMNGKRRLGFSFNHGSMANAMPQALGAQATEPERQ VVAMCGDGGFMSMLMGDFLSVVQMKLPVKIVVFNNSVLGFVAMEMKAGGYLT DGTELHDTNFARIAEACGITGIRVEKASEVDEALQRAFSIDGPVLVDVVVAKEEL AIPPQIKLEQAKGFSLYMLRAIISGRGDEVIELAKTNWLR
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	poxB
<b>Expression Region</b>	1-572aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	66.0kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.