



Recombinant Human KeRatin, type I cytoskeletal 18 (KRT18), partial

Product Code	CSB-BP012519HU2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05783
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ETENRRLESKIREHLEKKGPQVRDWSHYFKIIEDLRAQIFANTVDNARIVLQIDN ARLAADDFRVKYETELAMRQSVENDIHGLRKKVIDDTNITRLQLETEIEALKEELL FMKKNHEEEVKGLQAQIASSGLTVEVDAPKSQDLAKIMADIRAQYDELARKNR EELDKYWSQQIEESTTVTTQSAEVGAAETTLTELRRTVQSLEIDLDSMRNLKA SLENSLREVEARYALQMEQLNGILLHLESELAQTRAEGQRQAQEYEALLNIKVK LEAEIATYRRLLEDGEDFNLDAL
Research Area	Cell Cycle
Source	Baculovirus
Target Names	KRT18
Expression Region	102-396aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.