



Recombinant Ba (E.coli -k12) Flagellin protein

Product Code	CSB-YP300968ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04949
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	AQVINTNSLSLITQNNINKNQSALSSSIERSLSSGLRINSAKDDAAGQAIANRFTS NIKGLTQAARNANDGISVAQTTEGALSEINNQLRVRELTQATTGTNSESDLS SIQDEIKSRLDEIDRVSGQTQFNGVNVLAKNGSMKIQV GANDNQTITIDLKQIDA KTLGLDGFVSVKNNDTVTTAPVTAFGATTTNNIKLTGITLSTEAATDTGGTNPA SIEGVYTDNGNDYYAKITGGDNDGKYAVTVANDGTVTMATGATANATVTDA NTTKATTITSGGTPVQIDNTAGSATANLGAVSLVKLQDSKGNDDTYALKDTN GNLYAADVNETTGAVSVKTITYTDSSGAASSPTAVKLGDDGKTEVVDIDGKT YDSADLNGGNLQTGLTAGGEALTAVANGKTTDPLKALDDAIASVDKFRSSLGA VQNRLDSAVTNLNNNTTTLNSEAQSRIQDADYATEVSNMSKAQIIQQAGNSVLA KANQVPQQVLSLLQG
Research Area	Others
Source	Yeast
Target Names	fliC
Expression Region	2-498aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.