



Recombinant PI (Arabidopsis thaliana (Mouse-ear cress)) L-ascorbate peroxidase 1,

Product Code	CSB-EP134794PI-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42738
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	KNYPTVSEDYKKAVEKCRRLRGLIAEKNCAPIMVRLAWHSAGTFDCQSRTG GPFGTMRFDAEQAHGANSIHALRLLDPIREQFPTISFADFHLQAGVVAVEVT GGPDIPFHPGREDKPQPPPEGRLPDATKGC DHLRDVFAKQMGLSDKDIVALS GAHTLGRCHKDRSGFEGAWTSNPLIFDNSYFKELLSGEKEGLLQLVSDKALLD DPVFRPLVEKYAADEDAFFADYAEAHEMMLSELGF
Research Area	Cardiovascular
Source	E.coli
Target Names	CM1
Expression Region	3-247aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.