




# Recombinant Human Sorting nexin-20 (SNX20), partial

<b>Product Code</b>	CSB-EP022366HU
<b>Relevance</b>	Serves as a sorting protein that cycles P-selectin glycoprotein ligand 1 (PSLG1) into endosomes with no impact on leukocytes recruitment.
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z614
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MASPEHPGSPGCMGPITQCTARTQQEAPATGPDLPHPGPDGHLDTHTSGLSS NSSMTTRELQQYWQNQKCRWKHVKLLFEIASARIEERKVKVSKFVVYQIIVIQ
<b>Research Area</b>	Transport
<b>Source</b>	E.coli
<b>Target Names</b>	SNX20
<b>Expression Region</b>	1-102aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	38.4kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.