






Recombinant Pasteurella haemolytica Leukotoxin (lktA), partial

Product Code	CSB-EP314743ESE
Relevance	Pasteurella leukotoxins are exotoxins that attack host leukocytes and especially polymorphonuclear cells, by causing cell rupture. The leukotoxin binds to the host LFA-1 integrin and induces a signaling cascade leading to many biological effects, including tyrosine phosphorylation of the CD18 tail, elevation of the intracellular Ca ²⁺ and lysis of the host cell . This leukotoxin is a major contributor to the pathogenesis of lung injury in ovine pneumonic pasteurellosis. It has also weak holytic activity.
Abbreviation	Recombinant Mannheimia haemolytica lktA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C085
Product Type	Recombinant Protein
Immunogen Species	Mannheimia haemolytica (Pasteurella haemolytica)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	IGTSHNDIFKGSKFND AFNGGDGVDTIDGNDGNDRLFGGKGDDIIDGGNGDDF IDGGKGNLHGGKGDDIFVHRQGDGND SITESEGNKLSFSDSNLKD LTFEK VNHHLVITNTKQEKVTIQNW FREA EFAKTIQNYVATRDDKIEEII GQNGERITSK QVDELIEKGNGKIAQSELTKVVDNYQLLKYSRDASNLDKLISSASAFTSSNDS RNVLASPTSM LDP SLSSIQFARAA
Research Area	Others
Source	E.coli
Target Names	lktA
Expression Region	715-953aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	53.0kDa
Protein Length	Partial
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.