




Recombinant Human Glycosyltransferase 8 domain-containing protein 2 (GLT8D2)

Product Code	CSB-EP009536HU
Abbreviation	Recombinant Human GLT8D2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H1C3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	KKVHKGTVPKNDADDESETPEELEEEEIPVVICAAAGRMGATMAAINSIYSNTDANILFYVVGLRNTLTRIRKWIEHSKLREINFKIVEFNPMVLKGGKIRPDSSRPELLQP LNFVRFYLP LLIHQHEKVIYLD DDDVIVQGGDIQELYDTT LALGHAAAFSDDCDLPS AQDINRLVGLQNTYMGYLDYRKKAIKDLGISPSTCSFNPGVIVANMTEWKHQRITKQLEKWMQKNVEENLYSSSLGGGVATSPMLIVFHGKYSTINPLWHIRHLGW NPDARYSEHFLQEAKLLHWNGRHKPWDFPSVHNDLWESWFVPDPAGIFKLN HHS
Research Area	Metabolism
Source	E.coli
Target Names	GLT8D2
Expression Region	25-349aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	41.2kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.