



PRF Antibody, HRP conjugated

Product Code	CSB-PA15107B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q58NA1
Immunogen	Recombinant Toxoplasma gondii Profilin protein (2-163AA)
Raised In	Rabbit
Species Reactivity	Toxoplasma gondii
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Profilin, PRF
Immunogen Species	Toxoplasma gondii
Research Area	Others
Target Names	PRF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.