



FGF21 Antibody, Biotin conjugated

Product Code	CSB-PA06404D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSA1
Immunogen	Recombinant Human Fibroblast growth factor 21 protein (29-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 21 (FGF-21), FGF21
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGF21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.