



TFRC Antibody, HRP conjugated

Product Code	CSB-PA07219B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02786
Immunogen	Recombinant Human Transferrin receptor protein 1 protein (94-754AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transferrin receptor protein 1 (TR) (TfR) (TfR1) (Trfr) (T9) (p90) (CD antigen CD71) [Cleaved into: Transferrin receptor protein 1, serum form (sTfR)], TFRC
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TFRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.