



# L1CAM Antibody, HRP conjugated

Product Code	CSB-PA07439B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32004
Immunogen	Recombinant Human Neural cell adhesion molecule L1 protein (995-1119AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neural cell adhesion molecule L1 (N-CAM-L1) (NCAM-L1) (CD antigen CD171), L1CAM, CAML1 MIC5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	L1CAM