



TNF Antibody, Biotin conjugated

Product Code	CSB-PA023955DD01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23383
Immunogen	Recombinant Sheep Tumor necrosis factor protein (78-234AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor (Cachectin) (TNF-alpha) (Tumor necrosis factor ligand superfamily member 2) (TNF-a) [Cleaved into: Tumor necrosis factor, membrane form (N-terminal fragment) (NTF); Intracellular domain 1 (ICD1); Intracellular domain 2 (ICD2); C-domain 1; C-domain 2; Tumor necrosis factor, soluble form], TNF, TNFA TNFSF2
Immunogen Species	Ovis aries (Sheep)
Research Area	Others
Target Names	TNF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.