



RBP1 Antibody, Biotin conjugated

Product Code	CSB-PA019480YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09455
Immunogen	Recombinant Human Retinol-binding protein 1 protein (70-112AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinol-binding protein 1 (Cellular retinol-binding protein) (CRBP) (Cellular retinol-binding protein I) (CRBP-I), RBP1, CRBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.