



# fre Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA18537D0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0AEN1
<b>Immunogen</b>	Recombinant Escherichia coli NAD(P)H-flavin reductase protein (2-233AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	NAD(P)H-flavin reductase (EC 1.16.1.3) (EC 1.5.1.41) (Aquacobalamin reductase) (FMN reductase) (Ferrisiderophore reductase C) (NAD(P)H:flavin oxidoreductase) (Riboflavin reductase [NAD(P)H]), fre, fadI flrD fsrC ubiB
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Research Area</b>	Others
<b>Target Names</b>	fre
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.