




Recombinant Erythrina caffra Trypsin inhibitor DE-3

Product Code	CSB-YP362840EGO
Relevance	Inhibition of trypsin.
Abbreviation	Recombinant Erythrina caffra Trypsin inhibitor DE-3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09943
Product Type	Recombinant Protein
Immunogen Species	Erythrina caffra (Kaffir tree) (Coastal coral tree)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VLLDGNGEVQNGGTYYLLPQVWAQGGGVQLAKTGEETCPLTVVQSPNELS DGKPIRIESRLRSFIPDDDKVRIGFAYAPKCAPSPWWTVVEDEQEGLSVKLSE DESTQFDYPFKFEQVSDQLHSYKLLYCEGKHEKCASIGINRDQKGYRRLVTE DYPLTVVLKKDESS
Research Area	Others
Source	Yeast
Protein Names	Recommended name: Trypsin inhibitor DE-3
Expression Region	1-172aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.3kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.