



# blaNDM-1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA07617D0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	C7C422
<b>Immunogen</b>	Recombinant Klebsiella pneumoniae Metallo-beta-lactamase type 2 protein (29-270AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Klebsiella pneumoniae
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Metallo-beta-lactamase type 2 (EC 3.5.2.6) (B2 metallo-beta-lactamase) (Beta-lactamase type II) (Metallo-beta-lactamase NDM-1) (Metallo-beta-lactamase type II) (New Delhi metallo-beta-lactamase-1) (NDM-1), blaNDM-1
<b>Immunogen Species</b>	Klebsiella pneumoniae
<b>Research Area</b>	Others
<b>Target Names</b>	blaNDM-1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.