



DLL3 Antibody, Biotin conjugated

Product Code	CSB-PA882142LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYJ7
Immunogen	Recombinant Human Delta-like protein 3 protein (27-492AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delta-like protein 3 (Drosophila Delta homolog 3) (Delta3), DLL3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DLL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.