



Lgals7 Antibody

Product Code	CSB-PA09877A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54974
Immunogen	Recombinant Mouse Galectin-7 protein (2-136AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-7 (Gal-7), Lgals7
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Lgals7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.