



FCGRT Antibody, HRP conjugated

Product Code	CSB-PA851443LB01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8SPV9
Immunogen	Recombinant Macaca fascicularis IgG receptor FcRn large subunit p51 protein (24-297AA)
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IgG receptor FcRn large subunit p51 (FcRn) (IgG Fc fragment receptor transporter alpha chain) (Neonatal Fc receptor), FCGRT, FCRN
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Research Area	Others
Target Names	FCGRT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.