



GLP1R Antibody, HRP conjugated

Product Code	CSB-PA009514YB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43220
Immunogen	Recombinant Human Glucagon-like peptide 1 receptor protein (24-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucagon-like peptide 1 receptor (GLP-1 receptor) (GLP-1-R) (GLP-1R), GLP1R
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GLP1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.