



rodA Antibody, FITC conjugated

Product Code	CSB-PA340427LC01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41746
Immunogen	Recombinant Neosartorya fumigata Hydrophobin protein (19-159AA)
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydrophobin (Rodlet protein), rodA, hyp1
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Research Area	Others
Target Names	rodA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.