



GRA1 Antibody

Product Code	CSB-PA319538YA01TOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13403
Immunogen	Recombinant Toxoplasma gondii Dense granule protein 1 protein (26-187AA)
Raised In	Rabbit
Species Reactivity	Toxoplasma gondii
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Toxoplasma gondii
Research Area	Others
Target Names	GRA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.