



GRA1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA319538YC01TOV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P13403 |
| Immunogen | Recombinant Toxoplasma gondii Dense granule protein 1 protein (26-187AA) |
| Raised In | Rabbit |
| Species Reactivity | Toxoplasma gondii |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Dense granule protein 1 (Protein GRA 1) (Major antigen p24), GRA1 |
| Immunogen Species | Toxoplasma gondii |
| Research Area | Others |
| Target Names | GRA1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |