



GRA1 Antibody, Biotin conjugated

Product Code	CSB-PA319538YD01TOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13403
Immunogen	Recombinant Toxoplasma gondii Dense granule protein 1 protein (26-187AA)
Raised In	Rabbit
Species Reactivity	Toxoplasma gondii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dense granule protein 1 (Protein GRA 1) (Major antigen p24), GRA1
Immunogen Species	Toxoplasma gondii
Research Area	Others
Target Names	GRA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.