



CHI3L2 Antibody, Biotin conjugated

Product Code	CSB-PA620893LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15782
Immunogen	Recombinant Human Chitinase-3-like protein 2 protein (27-390AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chitinase-3-like protein 2 (Chondrocyte protein 39) (YKL-39), CHI3L2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	CHI3L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.