



Muc5ac Antibody

| | |
|---------------------|---|
| Product Code | CSB-PA17229A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | E9PWB6 |
| Immunogen | Recombinant Mouse Mucin 5, subtypes A and C, tracheobronchial/gastric protein (2452-2721AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mucin 5, subtypes A and C, tracheobronchial/gastric, Muc5ac |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Others |
| Target Names | Muc5ac |