



Muc5ac Antibody, Biotin conjugated

Product Code	CSB-PA17229D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E9PWB6
Immunogen	Recombinant Mouse Mucin 5, subtypes A and C, tracheobronchial/gastric protein (2452-2721AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mucin 5, subtypes A and C, tracheobronchial/gastric, Muc5ac
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Muc5ac
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.