



NEU3 Antibody, FITC conjugated

Product Code	CSB-PA17267C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ49
Immunogen	Recombinant Human Sialidase-3 protein (2-428AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialidase-3 (EC 3.2.1.18) (Ganglioside sialidase) (Membrane sialidase) (N-acetyl-alpha-neuraminidase 3), NEU3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NEU3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.