



Fgf23 Antibody, Biotin conjugated

Product Code	CSB-PA008629LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VI82
Immunogen	Recombinant Rat Fibroblast growth factor 23 protein (25-251AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 23 (FGF-23), Fgf23
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Fgf23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.