



Fgf23 Antibody, Biotin conjugated

Product Code	CSB-PA008629LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VI82
Immunogen	Recombinant Rat Fibroblast growth factor 23 protein (25-251AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Relevance	Regulator of phosphate homeostasis. Inhibits renal tubular phosphate transport by reducing SLC34A1 levels. Regulator of vitamin-D metabolism. Negatively regulates osteoblasts differentiation and matrix mineralization. Acts directly on the parathyroid to decrease PTH secretion. Upregulates EGR1 expression in the presence of KL.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 23 (FGF-23), Fgf23
Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Fgf23