



Dio3 Antibody

Product Code	CSB-PA344617LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49897
Immunogen	Recombinant Rat Thyroxine 5-deiodinase protein (37-278AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroxine 5-deiodinase (EC 1.21.99.3) (5DIII) (DIOIII) (Type 3 DI) (Type III iodothyronine deiodinase), Dio3, ltdi3 Txdi3
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Dio3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.