



IFIT3 Antibody, Biotin conjugated

Product Code	CSB-PA011022HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14879
Immunogen	Recombinant Human Interferon-induced protein with tetratricopeptide repeats 3 protein (1-490AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced protein with tetratricopeptide repeats 3 (IFIT-3) (CIG49) (ISG-60) (Interferon-induced 60 kDa protein) (IFI-60K) (Interferon-induced protein with tetratricopeptide repeats 4) (IFIT-4) (Retinoic acid-induced gene G protein) (P60) (RIG-G), IFIT3, CIG-49 IFI60 IFIT4 ISG60
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFIT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.