




# NABP1 Antibody

<b>Product Code</b>	CSB-PA850257HA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96AH0
<b>Immunogen</b>	Recombinant Human SOSS complex subunit B2 protein (1-204AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	SOSS complex subunit B2 (Nucleic acid-binding protein 1) (Oligonucleotide/oligosaccharide-binding fold-containing protein 2A) (Sensor of single-strand DNA complex subunit B2) (Sensor of ssDNA subunit B2) (SOSS-B2) (Single-stranded DNA-binding protein 2) (hSSB2), NABP1, OBFC2A SSB2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	NABP1

## Image



Western blot  
All lanes: NABP1 antibody at 2µg/ml  
Lane 1: Hela whole cell lysate  
Lane 2: Jurkat whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 23, 14, 15 kDa  
Observed band size: 23 kDa