



# USH1C Antibody

<b>Product Code</b>	CSB-PA896761HA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y6N9
<b>Immunogen</b>	Recombinant Human Harmonin protein (1-533AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Harmonin (Antigen NY-CO-38/NY-CO-37) (Autoimmune enteropathy-related antigen AIE-75) (Protein PDZ-73) (Renal carcinoma antigen NY-REN-3) (Usher syndrome type-1C protein), USH1C, AIE75
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	USH1C
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.