



Bag3 Antibody, HRP conjugated

Product Code	CSB-PA864372HB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JLV1
Immunogen	Recombinant Mouse BAG family molecular chaperone regulator 3 protein (2-577AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BAG family molecular chaperone regulator 3 (BAG-3) (Bcl-2-associated athanogene 3) (Bcl-2-binding protein Bis), Bag3, Bis
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Bag3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.