



PISD Antibody

Product Code	CSB-PA880077LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UG56
Immunogen	Recombinant Human Phosphatidylserine decarboxylase proenzyme protein (1-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylserine decarboxylase proenzyme, mitochondrial (EC 4.1.1.65) [Cleaved into: Phosphatidylserine decarboxylase beta chain; Phosphatidylserine decarboxylase alpha chain], PISD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PISD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.