



Histidine-rich protein PFHRP-II Monoclonal Antibody

Product Code	CSB-MA3564761A0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05227
Immunogen	Recombinant Histidine-rich protein PFHRP-II for Plasmodium falciparum
Raised In	mouse
Species Reactivity	Plasmodium falciparum
Specificity	specific for PFHRP-II denatured form
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	PFHRP-II
Product Type	Monoclonal Antibody
Immunogen Species	Plasmodium falciparum
Clone No.	2A1B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.